

# **SFWRENG 3K04: Software Development**

## **Assignment 2**

### **Simulink**

Group 25

December 1st, 2021

## 2. Simulink

### 2.1 Simulink Overview

### 2.2 Parameters

#### 2.2.1 Programmable Input Parameters Regarding Chart (Inputs)

### 2.3 Design Decision

#### 2.3.1 Overview

#### 2.3.2 AOO and VOO Design

#### 2.3.3 AAI and VVI Design

#### 2.3.4 DOO Design

#### 2.3.5 AOOR, VOOR, AAIR, VVIR, DOOR Design

#### 2.3.6 Subsystem Usage

### 2.4 State Diagram

#### 2.4.1 Input Chart

#### 2.4.2 Adaptive Rate Adjusting Chart

#### 2.4.3 Main Pacemaker Chart

### 2.5 Function Blocks

#### 2.5.1 Overview of the Accelerometer Data Processing Function Group

#### 2.5.2 Function Block: Convert vectors into magnitude

#### 2.5.3 Function Block: Filtering Coefficient

#### 2.5.4 Function Block: Calculate Target Rate

#### 2.5.5 Function: send\_parameters()

### 2.6 Hardware Hiding

#### 2.6.1 Input\_subsystem

#### 2.6.2 Pulse Detect subsystem

#### 2.6.3 Output\_subsystem

### 2.7 Testing

#### 2.7.1 Test case 1

#### 2.7.2 Test case 2

#### 2.7.3 Test case 3

#### 2.7.4 Test case 4

#### 2.7.5 Test case 5

### 2.8 Pin Map

### 2.9 Requirements Changes That are Likely

### 3.0 Design Decision That Are Likely To Change

## 2. Simulink

### 2.1 Simulink Overview

The purpose of the simulink implementation is to simulate the functionalities of a pacemaker.

The pacemaker design consists of 10 working states in total, which are VOO, AOO, VVI, VOO, VOOR, AOOR, VVIR, AAIR, DOO, and DOOR. VOO, AOO and DOO are states that will apply direct pacing, VVI and AAIR are states that need to sense the heart activities before determining whether to pace or not. VOOR, AOOR, VVIR, AAIR are modes that require adaptive rate pacing, which means the pacing rate depends on the activity level sensed by the accelerometer on the pacemaker.

### 2.2 Parameters

#### 2.2.1 Programmable Input Parameters Regarding Chart (Inputs)

| Name      | Initial Values | Units | Data Type | Functionality   |
|-----------|----------------|-------|-----------|---|
| Mode      | 1              |       | uint8     | Select from mode 1 - 10   |
| LRL       | 60             | ppm   | uint8     | Default pacing pulse rate   |
| MSR       | 160            | ppm   | uint8     | Maximum pacing rate allowed for adaptive pacing modes                           |
| ARP       | 250            | msec  | uint16    | Atrial Refractory Period  |
| VRP       | 250            | msec  | uint16    | Ventricle Refractory Period   |
| ATR_AMP   | 3.3            | Volt  | single    | Desired pacing PWM for atrial, the value will be convert to duty cycle later    |
| VENT_AMP  | 3.3            | Volt  | single    | Desired pacing PWM for ventricle, the value will be convert to duty cycle later |
| ATR_Width | 10             | msec  | uint8     | Pacing width for atrial pacing  |

|                  |     |      |        |  |
|------------------|-----|------|--------|--|
| VENT_Width       | 10  | msec | uint8  | Pacing width for ventricular pacing  |
| ATR_Sensitivity  | 3   | Volt | single | Threshold voltage for sensing a natural atrial beat.   |
| VENT_Sensitivity | 3   | Volt | single | Threshold voltage for sensing a natural ventricular beat.  |
| Threshold        | 1.1 |      | single | Threshold value for adaptive pacing mode. It decides when the LRL should increase. Range 1 - 6.9 |
| RF               | 16  |      | uint8  | Response factor, which decides the desired pacing rate.  |
| Reaction_Time    | 10  | sec  | uint8  | The time takes to increase from LRL to the MSR.  |
| Recovery_Time    | 30  | sec  | uint8  | The time takes to decrease from MSR to LRL.  |
| AV_delay         | 150 | msec | uint16 | AV delays for dual chamber pacing modes  |
| URL              | 160 | ppm  | uint8  | Upper rate limit for pacing  |

## 2.3 Design Decision

### 2.3.1 Overview

The simulink design consists of 6 parts in general. The input chart, the input subsystem, the main pacemaker, the data processing MatLab function group, the output subsystem and the set\_parameter() function block. By dividing the design into smaller parts, the dataflow relationship is more clear. The input data is being set in the input chart, then it flows into two directions, one to the set\_parameter() function block where the data is sent back to the DCM, the other is to the input subsystem where some data pre-processing occurs. From the input subsystem, while part of the data flows into the MatLab function group to perform some calculations, the other part together with the calculation results from the function group flow into the main pacemaker chart. Finally, the data flows to the output subsystem.

### 2.3.2 AOO and VOO Design

In these two modes, the design will continue pacing the atrial/ventricle. It is achieved by switching between the Charging\_Discharging state and the Pacing state. The Charging\_Discharging state integrates the functionalities to both charge/discharge the capacitor. The Pacing state simply paces the atrial/ventricle.

### 2.3.3 AAI and VVI Design

Firstly, the design will first charge C22 / discharge the blocking capacitor. Then, a state named Sensing\_VENT/ATR is inserted to determine if the ATR/VENT\_CMP\_DETECT detects a valid heartbeat. If the heart activity is sensed, it will loop back to the initial Charging\_Discharging state to prepare for the next detection, otherwise, it will proceed to the pacing state to pace the Atrial / Ventricle. After that, it will loop back to the initial Charging\_Discharging state.

### 2.3.4 DOO Design

The DOO mode consists of 5 states in total. The states are, Charging\_Discharging\_A, Pacing\_Atrial, ChargingV\_DischargingA, Pacing\_Ventricle, ChargingA\_DischargingV. The first Charging\_Discharging\_A state is a lead in state where it will only be executed one time after entry, then the code will loop among the other 4 states in the order that is mentioned above. It will allow the pacing to occur in the desired sequence with proper AV delay.

### 2.3.5 AOOR, VOOR, AAIR, VVIR, DOOR Design


The adaptive rate modes share the same design logic as their corresponding normal modes mentioned above. The only difference is that the pacing period in the rate adaptive mode depends on the calculated rate rather than the LRL.

### 2.3.6 Subsystem Usage

The main purpose of the subsystem is to implement the hardware hiding, the input/output pin assignments are all done in the subsystem. Also, the subsystem also makes the design look more clean when some basic data processing is required, such as data type conversion, all the data are calculated based on double type.


## 2.4 State Diagram

### 2.4.1 Input Chart



The Initialize state is used to set the initial value for all the parameters. And it will enter the standby state to determine whether to enter SET\_PARAM state or the ECHO\_PARAMETERS state. SET\_PARAM state is used to receive the serial input and assign the data to the corresponding variables. ECHO\_PARAMETERS is used to transfer all the data back to the DCM side by calling the function send\_parameters(). (the function will be introduced in the ‘2.5 Function blocks’ section)

## 2.4.2 Adaptive Rate Adjusting Chart



This state takes LRL, MSR, Reaction time, Recovery time, Target\_rate and Target\_LVL as input, variable ‘Target\_LVL’ is the activity level determined by the previous MatLab function block,


variable 'Target\_rate' is the target rate of the adaptive mode which is also determined by the previous MatLab function block. These function blocks will be introduced in the '2.5 Function Blocks' section later in this report. After processing, the state will return an output represented by variable 'Rate', which is the pacing rate used in the adaptive pacing mode. This rate is a dynamic variable, it will keep changing as the accelerometer data change.

The CMP\_LVL is an internal parameter, it represents the current activity level of the pacemaker, its value will keep updating as the activity level changes. After its initial value is set, it represents the current activity level and is used to compare with the incoming activity level to decide if the pacing rate needs to increase or decrease, then the CMP\_LVL is updated with that value and wait for the next comparison.

The chart starts with the entry state Entry\_only, where Rate and CMP\_LVL are set to default value (CMP\_LVL = 0, Rate = LRL). Then it enters a Stand\_by state, the code will compare the CMP\_LVL with the Target\_LVL. If (CMP\_LVL < Target\_LVL), it will enter the INCREASE state. If (CMP\_LVL > Target\_LVL), it will enter the DECREASE state. If (CMP\_LVL == Target\_LVL), it will wait at the Stand\_by state. In INCREASE state, the value of Rate will increment with the formula every 1 sec:  $\text{Rate} = \text{Rate} + ((\text{MSR} - \text{LRL}) / \text{Reaction\_time})$ , where '(MSR - LRL)/Reaction\_time' is the increment step, the increment will stop when the Rate reach Target\_rate. As long as CMP\_LVL is not equal to Target\_LVL, the state will transfer back to Stand\_by state for further decision. In DECREASE state, the logic is similar, the only difference is the formula,  $\text{Rate} = \text{Rate} - ((\text{MSR} - \text{LRL}) / \text{Recovery\_time})$ , the Rate will decrease until it reaches the Target\_rate. As long as CMP\_LVL is not equal to Target\_LVL, the state will transfer back to Stand\_by state for further decision.

### 2.4.3 Main Pacemaker Chart


#### (1). Overview




The entire process starts from the Entry\_Point, where the desired pacing mode is determined by the value of input mode. There are 10 modes that could be selected.

## (2). AOO and VOO. AOOR and VOOR

### AOO mode



### AOOR mode



AOO and VOO modes share similar logic. On the LHS, the Charging\_Discharging state integrates the function of charging for the next pace and discharge after pacing occurs. On the RHS, the Pacing state executes the pacing after the capacitor has been charged. In these two modes, the pacing period are both determined by LRL.



AOOR and VOOR have exactly the same design diagrams as their OO modes. The only difference is the pacing period, the pacing period of the R mode is determined by the calculated Adaptive rate (represented by variable: Rate) rather than LRL.

### (3). AAI and VVI. AAIR and VVIR

#### AAI mode



#### AAIR mode




AAI and VVI modes share similar logic. On the LHS, the Charging\_Discharging state integrates the function of charging for the next pace and discharge after pacing occurs. The Sensing\_ATR/VENT state creates the pacing inhibited time interval after natural beats have been sensed (It is when ATR/VENT\_CMP\_DETECT are true). On the RHS, the Pacing state executes the pacing after the capacitor has been charged and the inhibit time period has passed. In these two modes, the pacing period are both determined by LRL.


AAIR and VVIR have exactly the same design diagrams as their 'I' modes. The only difference is the pacing period, the pacing period of the R mode is determined by the calculated Adaptive rate (represented by variable: Rate) rather than LRL.

### (4) DOO and DOOR

#### DOO mode



## DOOR mode




The DOO mode consists of 5 states in total. However, the entry state Charging\_Dischaging\_A is a boost up state which will only be executed one time after entering this mode, the code will loop among the other 4 states. The Charging\_Dischaging\_A state charges the capacitor for Atrial pacing, after(pacing\_period - atr\_width) msec, the Pacing\_Atrial state executes the atrial pacing, then after (atr\_width) msec, ChargingV\_DischargingA state discharges for atrial and recharges for the next ventricular pacing. After(av\_delay - atr\_width) msec, Pacing\_Ventricle state paces the ventricle, then after(vent\_width) msec, the ChargingA\_DischargingV state discharges for the ventricle and recharges the capacitor for next atrial pacing, after(pacing\_period - av\_delay - vent\_width) msec, the Pacing\_Atrial state paces the atrial again. For the rest of the time, the code will loop among the above 4 states. For this DOO mode, all the pacing period depends on the LRL.

The DOOR mode has the same design logic as the DOO mode, the only difference is the pacing period depends on the calculated adaptive pacing rate (represented by variable 'Rate').

## 2.5 Function Blocks

### 2.5.1 Overview of the Accelerometer Data Processing Function Group



This function group takes the raw data from the accelerometer and outputs the calculated activity level and the corresponding target rate.

### 2.5.2 Function Block: Convert vectors into magnitude

```

1  %This function calculate the magnitude of the input-3dimentional-vector input
2  function Processed_data = fcn(Raw_data)
3
4      Processed_data = sqrt((Raw_data(1))^2+(Raw_data(2))^2+(Raw_data(3))^2);

```

This function calculates the magnitude of the accelerometer data, since the data represents 3 values from the x, y, z axis.

### 2.5.3 Function Block: Filtering Coefficient

```

1  %This function creates the coefficient for the signal filter.
2  %It generates a array of all ones, it takes each 200 samples and calculate
3  %the average.
4  function coefficient = fcn()
5      coefficient = (1/200)*ones(1,200);
6

```

This function creates the coefficient for the signal filter. It generates an array of all ones. The idea is to take each 200 samples and calculate the average.

### 2.5.4 Function Block: Calculate Target Rate

```

1  %This function calculates the target pacing rate corresponding to the activity
2  %level.
3  function Target_rate = fcn(LRL, Response_factor, MSR, Activity_level, Threshold)
4
5      if(Activity_level > Threshold)
6
7          x1 = LRL + Activity_level*Response_factor;
8          if(x1 > MSR)
9              Target_rate = MSR;
10         else
11             Target_rate = x1;
12         end
13     else
14         Target_rate = LRL;
15     end


```

The Activity\_level is the output value from the FIR filter. The function decides the target rate by adding the multiplication result of the activity level and the response factor to the base LRL.

When the calculated value exceeds MSR, then it is clipped to the MSR value.

The discussed calculation will only execute if the activity level exceeds the input activity level, otherwise, the target rate will equal to the LRL.


### 2.5.5 Function: send\_parameters()



When this function is called in the Input state chart, it will send all the parameters back to the DCM.


## 2.6 Hardware Hiding

### 2.6.1 Input\_subsystem




This subsystem isolates the input chart from the implementation. In this intermediate subsystem, the input parameters can be adjusted according to the actual implementation needs. In this pacemaker implantation, all the calculations are done in type Double, and the input voltage is converted to duty cycle here.

## 2.6.2 Pulse Detect subsystem



This subsystem contains the input pins for ATR\_CMP\_DETECT and VENT\_CMP\_DETECT.

## 2.6.3 Output\_subsystem



This output subsystem groups all the output pins.

## 2.7 Testing

| Test Case                                   | Modes select            | Input Parameters   | Expected output                           | Actual Output                             | Pass / Fail |
|---|-------------------------|--|---|---|-------------|
| 1. Shake the pacemaker from still condition | DOOR,<br>AOOR,<br>VOOR, | LRL = 60;<br>MSR = 160<br>ARP = 250<br>VRP = 250<br>ATR_Amp = 3.3<br>VENT_Amp = 3.3<br>ATR_Width = 10<br>VENT_Width = 10<br>ATR_Sensitivity = 3<br>VENT_Sensitivity = 3<br>Threshold = 1.1<br>RF = 16<br>Reaction_time = 10<br>Recovery_time = 30<br>AV_delay = 150<br>URL = 160 | Pacing rate increase until it reaches MSR | Pacing rate increase until it reaches MSR | Pass        |
| 2. Rest the pacemaker after shaking         | DOOR,<br>AOOR,<br>VOOR  | LRL = 60;<br>MSR = 160<br>ARP = 250<br>VRP = 250<br>ATR_Amp = 3.3<br>VENT_Amp = 3.3<br>ATR_Width = 10<br>VENT_Width = 10<br>ATR_Sensitivity  | Pacing rate decrease until it reaches LRL | Pacing rate decrease until it reaches LRL | Pass        |


|  |            |  |  |  |      |
|--|------------|--|--|--|------|
|  |            | ity = 3<br>VENT_Sensitivity = 3<br>Threshold = 1.1<br>RF = 16<br>Reaction_time = 10<br>Recovery_time = 30<br>AV_delay = 150<br>URL = 160   |  |  |      |
| 3. Dual chamber pacing modes react correctly to the AV delay | DOO/DOOR   | LRL = 30;<br>MSR = 175<br>ARP = 250<br>VRP = 250<br>ATR_Amp = 3.3<br>VENT_Amp = 3.3<br>ATR_Width = 10<br>VENT_Width = 10<br>ATR_Sensitivity = 3<br>VENT_Sensitivity = 3<br>Threshold = 1.1<br>RF = 16<br>Reaction_time = 10<br>Recovery_time = 30<br>AV_delay = 200<br>URL = 175 | There is correct AV delay interval between atrial and ventricular pacing | There is correct AV delay interval between atrial and ventricular pacing | Pass |
| 4. Shake the pacemaker from still condition                  | AAIR, VVIR | LRL = 60;<br>MSR = 170<br>ARP = 240<br>VRP = 250<br>ATR_Amp =  | Pacing rate increases while still maintain the required                  | Pacing rate increases while still maintain the required                  | Pass |



|                                     |            |   |   |   |      |
|-------------------------------------|------------|---|---|---|------|
|                                     |            | 3.3<br>VENT_Amp = 3.3<br>ATR_Width = 10<br>VENT_Width = 10<br>ATR_Sensitivity = 3<br>VENT_Sensitivity = 3<br>Threshold = 1.1<br>RF = 12<br>Reaction_time = 15<br>Recovery_time = 40<br>AV_delay = 150<br>URL = 170  | inhibition time interval after natural beat is sensed   | inhibition time interval after natural beat is sensed   |      |
| 5. Rest the pacemaker after shaking | AAIR, VVIR | LRL = 60;<br>MSR = 165<br>ARP = 250<br>VRP = 250<br>ATR_Amp = 3.3<br>VENT_Amp = 3.3<br>ATR_Width = 10<br>VENT_Width = 10<br>ATR_Sensitivity = 3<br>VENT_Sensitivity = 3<br>Threshold = 1.1<br>RF = 16<br>Reaction_time = 16<br>Recovery_time = 60<br>AV_delay = | Pacing rate decreases while still maintain the required inhibition time interval after natural beat is sensed | Pacing rate decreases while still maintain the required inhibition time interval after natural beat is sensed | Pass |


|  |  |                  |  |  |  |
|--|--|------------------|--|--|--|
|  |  | 150<br>URL = 165 |  |  |  |
|--|--|------------------|--|--|--|

2.7.1 Test case 1




The example pictures were taken under DOOR mode. The first picture shows that the pacemaker is resting on the table, the second one shows that when the pacemaker is being shaken. It is obvious that the pacing rate is getting faster when the pacemaker is on the move.

2.7.2 Test case 2






The example pictures were taken under DOOR mode. Following test case 1, when the pacemaker is resting on the table after going through activity, it is obvious that the pacing rate is slower.

### 2.7.3 Test case 3




The example picture was taken under DOO mode. The AV delay is set to be 200 msec, it is obvious that the order of pacing and time intervals between each pacing are correct, according to the input settings.

## 2.7.4 Test case 4



The example pictures were taken under AAIR mode. From the first two pictures, it can tell that when the pacemaker is resting on the table, it shows the same habitat as the normal sensing mode. The third picture shows that when the pacemaker is being shaken, the pacing rate increases while still maintaining the required inhibition time interval after natural beat is sensed.

## 2.7.5 Test case 5



The example pictures were taken under VVIR mode. The first picture is when the pacemaker is being shaken, the second one is when the pacemaker comes to a rest. It is obvious that the pacing rate is slower, and at the same time, it still maintains the required inhibition time interval after natural beat is sensed.

## 2.8 Pin Map

| Input            |           |                                 |  |
|------------------|-----------|---------------------------------|--|
| Name             | Pin_Name  | Connection                      | Functionality  |
| ATR_CMP_DETECT   | PTC16(D0) | ATR_Detect                      | Used To check if the signal voltage is higher than threshold voltage which will make it output ON to start the atrium sensing. |
| VENT_CMP_DETECT  | PTC17(D1) | VENT_Detect                     | Same functionality as in ATR_CMP_DETECT but for the ventricle sensing signal   |
| ATR_SIGNAL       | PTB2 (A0) | FRDM-K64F Serial Transmit Block | Send the atrial signal to the DCM  |
| VENT_SIGNAL      | PTB3 (A1) | FRDM-K64F Serial Transmit Block | Send the ventricle signal to the DCM   |
| Output           |           |                                 |  |
| Name             | Pin_Name  | Connection                      | Functionality  |
| PACE_CHARGE_CTRL | PTB9(D2)  | pace_charge_ctrl                | To charge C22 when it is ON, and discharge C22 when it is off  |
| VENT_CMP_REF_PWM | PTA1(D3)  | vent_cmp_ref_pwm                | To set a threshold voltage when the ventricle is going to be sensed  |

|                 |           |                 |  |
|-----------------|-----------|-----------------|--|
| Z_ATR_CTRL      | PTB23(D4) | z_atr_ctrl      | Allows the impedance circuit to be connected to the ring electrode of the atrium |
| PACING_REF_PWM  | PTA2(D5)  | pacing_ref_pwm  | Used to charge C22   |
| ATR_CMP_REF_PWM | PTC2(D6)  | atr_cmp_ref_pwm | Same functionality as in VENT CMP REF PWM but for the atrial action potential.   |
| Z_VENT_CTRL     | PTC3(D7)  | z_vent_ctrl     | Used identically to Z ATR CTRL but for the ventricle.                            |
| ATR_PACE_CTRL   | PTA0(D8)  | atr_pace_ctrl   | Discharge C22 through the atrium.  |
| VENT_PACE_CTRL  | PTC4(D9)  | vent_pace_ctrl  | Same functionality as in ATR PACE CTRL but for the ventricle.                    |
| PACE_GND_CTRL   | PTD0(D10) | pace_gnd_ctrl   | Used to connect controls the switch directly following the tip                   |
| ATR_GND_CTRL    | PTD2(D11) | atr_gnd_ctrl    | Used to connect the ATR RING OUT to GND.   |
| VENT_GND_CTRL   | PTD3(D12) | vent_gnd_ctrl   | Same functionality as in ATR_GND_CTRL but for the ventricle.                     |

|               |           |             |   |
|---------------|-----------|-------------|---|
| FRONTEND_CTRL | PTD1(D13) | sensingCTRL | Used to activate the sensing circuitry. The circuit outputs high when it is true. |
|---------------|-----------|-------------|---|

## 2.9 Requirements Changes That are Likely

1. Implement more modes.
2. More threshold values could be assigned if there are more modes

## 3.0 Design Decisions That Are Likely To Change

1. If more threshold values are assigned, the MatLab function group may need to be modified, since different calculating algorithms will be required.
2. More states will be added to the main pacemaker chart, if more modes are implemented.